



Test and Bottle

Test	Bottle
Microbiology	
Drinking Water (presence/absence)	120 mL Sterile
Fecal Coliform (MTF)	120 mL Sterile
Inorganic Chemistry	
Nitrate	250 mL
Arsenic	250 mL
Metals (Iron, Mg, Mn)	250mL to 1 liter Depending on number of metals
Mercury	250 mL
Whatcom Short List	1 Liter
Conductivity	250mL
pH	250mL
Complete Inorganics	2- 1 liter
Oil & Grease	HCL Brown Glass
TSS	250 mL
BOD	1 liter
Total Phosphorus	250 mL
Lead/Copper	Cube
Organic Chemistry	
VOC (Drinking Water)	Special vials (please call)
VOC (Wastewater)	Special vials (please call)
HAA/TTHM	Special vials (please call)
NWTPH-DX	Special vials (please call)
SOC Package (WA)	Special vials (please call)
Herbicides/Pesticides	Special vials (please call)